

# Silicon Analysis in Petroleum and Bio Fuels

From gasoline to ethanol and toluene, Signal<sup>®</sup> delivers total silicon analysis. Powered by MWDXRF<sup>®</sup>, Signal is a robust analysis solution for demanding petroleum and industrial environments.

#### **Applications**

- Total silicon analysis in hydrocarbons and bio fuels
- For use in refinery labs, pipeline terminals, additive plants, and inspection laboratories

#### **Features and Benefits**

- LOD: 0.65 ppm at 600 s
- Dynamic Range: 0.65 ppm 3000 ppm
- · Fits on any lab bench
- Touch Screen user interface
- User programmable measurement time: 30-900 s
- No conversion gasses, heating elements, quartz tubes or columns
- 75 W air-cooled excitation tube

#### Options

· LIMS data output compatible software





## **TRUSTED PRECISION**

Monochromatic Wavelength Dispersive X-ray Fluorescence (MWDXRF) utilizes state-of-the-art focusing and monochromating optics to increase excitation intensity and dramatically improve signalto-background over high power traditional WDXRF instruments. This enables significantly improved detection limits and precision, and a reduced sensitivity to matrix effects. A monochromatic and focused primary beam excites the sample and secondary characteristic fluorescence X-rays are emitted from the sample. A second monochromating optic selects the silicon characteristic X-rays and directs these X-rays to the detector. MWDXRF is a direct measurement technique and does not require consumable gasses or sample conversion.



## Low Range Calibration



## **Product Specifications**

-		
Model	Signal	
Test Method	ASTM D7757	
Dimensions	37 cm (w) x 50 cm (d) x 34 cm (h)	
Power	100-120 VAC, 47-63 HZ at 6.0 Amps/ 200-240 VAC, 47-63 HZ at 6.0 Amps	
Sample Cup Volume	10 ml	
Ambient Temperature Requirements	5-40° C (40-104° F)	
Dynamic Range	0.65 ppm - 3000 ppm	
Measurement	User selectable: 30-900 s	
Calibration	8 calibration curves. Automatic and manual calibration functionality	

### LOW RANGE CALIBRATION

Signal uses a weighted least squares regression which is extremely linear and easy to set up. Typical correlation (R value) is expected to be on the order of 0.999 or better.

**Precision** Typical repeatability (r) and reproducibility (R) values in gasoline, at 95% confidence. 600 s measurement time.

Silicon Concentration (ppm)	r	R
2	0.4	0.7
5	0.5	0.8
8	0.6	1.0
15	0.8	1.4
100	2	4
500	5	10



15 Tech Valley Drive, East Greenbush, New York 12061 USA PH 518.880.1500 FAX 518.880.1510 info@xos.com xos.com

